


MICROAIRE[®]

Pulse Lavage | Blades | Burs | K-Wires | Pins



To order call **800.648.4309** or email inquiry@microaire.com.




Single-Use Pulse Lavage System

- Battery-operated (included)
- For use in total-joint procedures
- Concurrent suction with soft cone splash-shield for optimal fluid containment, without damaging soft tissue
- Ergonomic shape to reduce muscle fatigue
- Complete single-use system
- For total knee and hip arthroplasty, and trauma/wound debridement


201-03-00 KNEE KIT

- Handpiece and battery
- Shower spray tip with splash shield



201-04-00 HIP KIT

- Handpiece and battery
- Shower spray tip with splash shield
- Femoral spray tip






201-05-00 FAN SPRAY KIT









- Handpiece and battery
- Fan spray tip with splash shield





Oscillating Saw Blades

Part No.	Width (mm)	Usable Length (mm)	Cut Thickness (mm)	Material Thickness (mm)	TPI	
ZO-060	5.8	12	0.6	0.4	36	
ZO-062	5.8	18.5	0.6	0.4	36	
ZO-063	9.7	18.5	0.6	0.4	36	

Sagittal Saw Blades







Part No.	Width (mm)	Usable Length (mm)	Cut Thickness (mm)	Material Thickness (mm)	TPI	
ZS-032	5.8	25.4	0.6	0.4	22	
ZS-033	14	41	0.6	0.4	22	
ZS-038	9.4	25.4	0.6	0.4	36	
ZS-039	4.5	25.4	0.6	0.4	36	
ZS-333	14	41	0.6	0.4	18	
ZS-338	9.4	25.4	0.5	0.4	34	
ZS-339	4.5	25.4	0.6	0.4	34	
ZS-343	9.5	25.4	0.6	0.4	18	

Sternum Saw Blades

Part No.	Width (mm)	Usable Length (mm)	Cut Thickness (mm)	Material Thickness (mm)	TPI	
ZR-032M	9.4	34	1.0	0.6	18	
ZR-033M	9.4	34	1.0	0.6	16	


Medium-Length Burs

Drill with Built-on Medium Bur Guard

Type	Part No.	Head Size Diameter (mm)	Head Length (mm)	# of Flutes	Overall Length (mm)	
Round-Cut Carbide	ZB-126	3	3	8	50.8	
	ZB-128	4	4	8	48.3	
	ZB-130	5	5	10	48.3	
Oval-Cut Carbide	ZB-136	4	7.9	8	47.8	
	ZB-138	5.5	9.9	10	47.8	
Wire Pass Drill	ZB-148	1.5	19.1	2	50.8	

Long Burs

All require long bur guard.

Type	Part No.	Head Size Diameter (mm)	Head Length (mm)	# of Flutes	Overall Length (mm)	
Oval-Cut Carbide	ZB-238	5.5	9.9	10	69.9	

Extra-Long Burs

All require extra bur guard.


Type	Part No.	Head Size Diameter (mm)	Head Length (mm)	# of Flutes	Overall Length (mm)	
Oval-Cut Carbide	ZB-336	4	7.9	8	95.3	



K-Wires

To order non-sterile, add "NS" to the end of the product number (e.g. 1600-XXXNS)

Length	Point Style	Type	Wire Diameter			
			0.9 mm / .035"	1.1 mm / .045"	1.4 mm / .054"	1.6 mm / .062"
4"	Single Trocar	Smooth		1604-145		1604-162
	Double Trocar	Smooth	1600-435	1600-445	1600-454	1600-462
6"	Double Trocar	Smooth	1600-635	1600-645	1600-654	1600-662
9"	Single Trocar	Smooth				1600-9625
	Double Trocar	Smooth		1600-945		1600-962




K-Wires with Wire Guide


Length	Point Style	Type	Wire Diameter	
			1.1 mm / .045"	1.6 mm / .062"
5.5"	Double Trocar	Smooth	1600-022	1600-023

K-Wire + Pin Protectors

Molded thermoplastic, protection from sharp points.

For K-Wires	Wire Diameter				Color	
	0.7 mm / .028"	0.9 mm / .035"	1.1 mm / .045"	1.6 mm / .062"		
	1601-028	1601-035	1601-045	1601-062	Yellow	
		Blue	White	Green		

For Steinmann Pins	Wire Diameter							Color
	2.0 mm / 5/64"	2.4 mm / 3/32"	2.8 mm / 7/64"	3.2 mm / 1/8"	3.6 mm / 9/64"	4.0 mm / 5/32"	4.8 mm / 3/16"	
	1601-564	1601-332	1601-764	1601-180	1601-964	1601-532	1601-316	Pink
		Gray	Tan	Lavendar	Brown	Navy Blue	Orange	



Steinmann Pins

To order non-sterile, add "NS" to the end of the product number (e.g. 1600-XXXNS)

Length	Point Style	Type	Wire Diameter					
			2.0 mm / 5/64"	2.4 mm / 3/32"	2.8 mm / 7/64"	3.2 mm / 1/8"	3.6 mm / 9/64"	4.8 mm / 3/16"
9"	Single Trocar	Smooth	1620-509	1624-509	1628-509	1632-509	1636-509	1648-509
	Double Trocar	Smooth	1620-109	1624-109				
	Single Diamond	Smooth	1620-609	1624-609				
	Double Diamond	Smooth	1620-709					